

For Research Use Only

AGXT Recombinant antibody, PBS Only (Detector)



Catalog Number: 83385-4-PBS

Basic Information

Catalog Number: 83385-4-PBS	GenBank Accession Number: BC132819	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 189	CloneNo.: 240305D2
Source: Rabbit	UNIPROT ID: P21549	
Isotype: IgG	Full Name: alanine-glyoxylate aminotransferase	
Immunogen Catalog Number: AG18062	Calculated MW: 392 aa, 43 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00420-2, AGXT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83385-5-PBS. Detection antibody: 83385-4-PBS. Standard: Ag18062. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00420-3, AGXT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83385-1-PBS. Detection antibody: 83385-4-PBS. Standard: Ag18062. Range: 1.25-80 ng/mL