

For Research Use Only

# AGXT Recombinant antibody, PBS Only (Capture)



Catalog Number: 83385-5-PBS

## Basic Information

Catalog Number: 83385-5-PBS	GenBank Accession Number: BC132819	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 189	CloneNo.: 240305F9
Source: Rabbit	UNIPROT ID: P21549	
Isotype: IgG	Full Name: alanine-glyoxylate aminotransferase	
Immunogen Catalog Number: AG18062	Calculated MW: 392 aa, 43 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00420-2, AGXT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83385-5-PBS. Detection antibody: 83385-4-PBS. Standard: Ag18062. Range: 0.625-80 ng/mL