

For Research Use Only

DSC3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83386-1-PBS

Basic Information

Catalog Number: 83386-1-PBS	GenBank Accession Number: NM_001941	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1825	CloneNo.: 240292A12
Source: Rabbit	UNIPROT ID: Q14574	
Isotype: IgG	Full Name: desmocollin 3	
Immunogen Catalog Number: AG34648	Calculated MW: 100 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00421-2, DSC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83386-2-PBS. Detection antibody: 83386-1-PBS. Standard: Ag34648. Range: 0.312-20 ng/mL