

For Research Use Only

# DSC3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83386-4-PBS

## Basic Information

Catalog Number: 83386-4-PBS	GenBank Accession Number: NM_001941	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1825	CloneNo.: 240292H4
Source: Rabbit	UNIPROT ID: Q14574	
Isotype: IgG	Full Name: desmocollin 3	
Immunogen Catalog Number: AG34648	Calculated MW: 100 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00421-1, DSC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83386-2-PBS. Detection antibody: 83386-4-PBS. Standard: Ag34648. Range: 0.312-40 ng/mL



Cytometric bead array standard curve of MP00421-3, DSC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83386-4-PBS. Detection antibody: 83386-3-PBS. Standard: Ag34648. Range: 0.625-20 ng/mL