

For Research Use Only

SIAH1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83389-1-PBS

Basic Information

Catalog Number: 83389-1-PBS	GenBank Accession Number: BC035562	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6477	CloneNo.: 240377A11
Source: Rabbit	UNIPROT ID: Q8IUQ4	
Isotype: IgG	Full Name: seven in absentia homolog 1 (Drosophila)	
Immunogen Catalog Number: AG34618	Calculated MW: 282 aa, 31 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00403-1, SIAH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83389-1-PBS. Detection antibody: 83389-3-PBS. Standard: Ag34618. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00403-2, SIAH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83389-2-PBS. Detection antibody: 83389-1-PBS. Standard: Ag34618. Range: 0.625-20 ng/mL