

For Research Use Only

SIAH1 Recombinant antibody, PBS Only



Catalog Number: 83389-4-PBS

Basic Information

Catalog Number: 83389-4-PBS	GenBank Accession Number: BC035562	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6477	CloneNo.: 240377B4
Source: Rabbit	UNIPROT ID: Q8IUQ4	
Isotype: IgG	Full Name: seven in absentia homolog 1 (Drosophila)	
Immunogen Catalog Number: AG34618	Calculated MW: 282 aa, 31 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83389-4-RR (SIAH1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83389-4-PBS in a different storage buffer formulation.



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 83389-4-RR (SIAH1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83389-4-PBS in a different storage buffer formulation.