

For Research Use Only

BRCA1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83390-3-PBS

Basic Information

Catalog Number: 83390-3-PBS	GenBank Accession Number: BC115037	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 672	CloneNo.: 240285D2
Source: Rabbit	UNIPROT ID: P38398	
Isotype: IgG	Full Name: breast cancer 1, early onset	
Immunogen Catalog Number: AG17946	Calculated MW: 1863 aa, 208 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00425-1, BRCA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83390-3-PBS. Detection antibody: 83390-5-PBS. Standard: Ag17946. Range: 0.313-40 ng/mL