

For Research Use Only

NOL11 Recombinant antibody, PBS Only (Capture)



Catalog Number:83391-3-PBS

Basic Information

Catalog Number: 83391-3-PBS	GenBank Accession Number: BC064404	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 25926	CloneNo.: 240371D2
Source: Rabbit	UNIPROT ID: Q9H8H0	
Isotype: IgG	Full Name: nucleolar protein 11	
Immunogen Catalog Number: AG11783	Calculated MW: 719 aa, 81 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00427-1, NOL11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83391-3-PBS. Detection antibody: 83391-2-PBS. Standard: Ag11783. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00427-3, NOL11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83391-3-PBS. Detection antibody: 83391-1-PBS. Standard: Ag11783. Range: 6.25-100 ng/mL