

For Research Use Only

Fas/CD95 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83400-4-PBS

Basic Information

Catalog Number:

83400-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000043.6

GeneID (NCBI):

355

UNIPROT ID:

P25445-1

Full Name:

Fas (TNF receptor superfamily,
member 6)

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

240340C5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00430-1, Fas/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83400-4-PBS. Detection antibody: 83400-1-PBS. Standard: Eg0877. Range: 0.625-40 ng/mL.