

For Research Use Only

FGFR3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83402-2-PBS

Basic Information

| | | |
|---------------------------------------|--|---|
| Catalog Number: 83402-2-PBS | GenBank Accession Number: GeneID (NCBI): 2261 | Purification Method: Protein A purification |
| Size: 1 mg/ml | UNIPROT ID: P22607-1 | CloneNo.: 240378G4 |
| Source: Rabbit | Full Name: fibroblast growth factor receptor 3 | |
| Isotype: IgG | Calculated MW: 88kd | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00407-1, FGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83402-1-PBS. Detection antibody: 83402-2-PBS. Standard: Eg1040. Range: 3.125-50 ng/mL