

For Research Use Only

# SLC2A6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83406-3-PBS

## Basic Information

<b>Catalog Number:</b> 83406-3-PBS	<b>GenBank Accession Number:</b> BC013740	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 11182	<b>CloneNo.:</b> 240388F4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UGQ3	
<b>Isotype:</b> IgG	<b>Full Name:</b> solute carrier family 2 (facilitated glucose transporter), member 6	
<b>Immunogen Catalog Number:</b> AG18481	<b>Calculated MW:</b> 55 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00414-1, SLC2A6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83406-1-PBS. Detection antibody: 83406-3-PBS. Standard: Ag18481. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00414-2, SLC2A6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83406-2-PBS. Detection antibody: 83406-3-PBS. Standard: Ag18481. Range: 3.13-100 ng/mL