

For Research Use Only

FEZ2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83407-5-PBS

Basic Information

Catalog Number: 83407-5-PBS	GenBank Accession Number: BC110386	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9637	CloneNo.: 240248G8
Source: Rabbit	UNIPROT ID: Q9UHY8	
Isotype: IgG	Full Name: fasciculation and elongation protein zeta 2 (zygin II)	
Immunogen Catalog Number: AG15370	Calculated MW: 353 aa, 40 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00422-1, FEZ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83407-2-PBS. Detection antibody: 83407-5-PBS. Standard: Ag15370. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00422-2, FEZ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83407-4-PBS. Detection antibody: 83407-5-PBS. Standard: Ag15370. Range: 0.313-40 ng/mL