

For Research Use Only

# Lactotransferrin/LTF Recombinant antibody, PBS Only (Capture)



Catalog Number: 83409-2-PBS

## Basic Information

Catalog Number:

83409-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002343.6

GeneID (NCBI):

4057

UNIPROT ID:

P02788-1

Full Name:

lactotransferrin

Calculated MW:

78 kDa

Purification Method:

Protein A purification

CloneNo.:

240280F9

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00429-2, Lactotransferrin/LTF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83409-2-PBS. Detection antibody: 83409-1-PBS. Standard: Eg0888. Range: 0.625-80 ng/mL