

For Research Use Only

C13orf18 Recombinant antibody, PBS Only (Capture)



Catalog Number:83411-2-PBS

Basic Information

Catalog Number:

83411-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15474

GenBank Accession Number:

BC032311

GeneID (NCBI):

80183

UNIPROT ID:

Q9H714

Full Name:

chromosome 13 open reading frame

18

Calculated MW:

662 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

240259G4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00398-2, C13orf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83411-2-PBS. Detection antibody: 83411-4-PBS. Standard: Ag15474. Range: 0.625-40 ng/mL.