

For Research Use Only

C13orf18 Recombinant antibody, PBS Only (Capture)



Catalog Number:83411-5-PBS

Basic Information

Catalog Number: 83411-5-PBS	GenBank Accession Number: BC032311	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 80183	CloneNo.: 240259F12
Source: Rabbit	UNIPROT ID: Q9H714	
Isotype: IgG	Full Name: chromosome 13 open reading frame 18	
Immunogen Catalog Number: AG15474	Calculated MW: 662 aa, 73 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00398-1, C13orf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83411-5-PBS. Detection antibody: 83411-3-PBS. Standard: Ag15474. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00398-3, C13orf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83411-5-PBS. Detection antibody: 83411-1-PBS. Standard: Ag15474. Range: 0.625-40 ng/mL.