

For Research Use Only

# Mouse Osteopontin Recombinant antibody, PBS Only (Capture)

Catalog Number: 83413-1-PBS



## Basic Information

<b>Catalog Number:</b> 83413-1-PBS	<b>GenBank Accession Number:</b> NM_001204233.1	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 20750	<b>CloneNo.:</b> 240250F7
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q547B5	
<b>Isotype:</b> IgG	<b>Full Name:</b> secreted phosphoprotein 1	
	<b>Calculated MW:</b> 32kd	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00400-1, MOUSE Osteopontin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83413-1-PBS. Detection antibody: 83413-3-PBS. Standard: Eg0317. Range: 0.313-5 ng/mL