

For Research Use Only

Cytokeratin 77 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83415-1-PBS



Basic Information

Catalog Number: 83415-1-PBS	GenBank Accession Number: BC122558	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 374454	CloneNo.: 240490A12
Source: Rabbit	UNIPROT ID: Q7Z794	
Isotype: IgG	Full Name: keratin 77	
Immunogen Catalog Number: AG21942	Calculated MW: 578 aa, 62 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00406-1, KRT77 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83415-2-PBS. Detection antibody: 83415-1-PBS. Standard: Ag21942. Range: 0.625-40 ng/mL