

For Research Use Only

Thrombomodulin Recombinant antibody, PBS Only (Detector)



Catalog Number:83416-2-PBS

Basic Information

Catalog Number:

83416-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000361.3

GeneID (NCBI):

7056

UNIPROT ID:

P07204

Full Name:

thrombomodulin

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

240365C1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00411-1, Thrombomodulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83416-1-PBS. Detection antibody: 83416-2-PBS. Standard: Eg0944. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00411-2, Thrombomodulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83416-3-PBS. Detection antibody: 83416-2-PBS. Standard: Eg0944. Range: 0.313-40 ng/mL