

For Research Use Only

# KIAA1715 Recombinant antibody, PBS Only (Capture)

Catalog Number:83418-2-PBS



## Basic Information

Catalog Number: 83418-2-PBS	GenBank Accession Number: BC105132	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 80856	CloneNo.: 240211G11
Source: Rabbit	UNIPROT ID: Q9C0E8	
Isotype: IgG	Full Name: KIAA1715	
Immunogen Catalog Number: AG34669		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00417-2, KIAA1715 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83418-2-PBS. Detection antibody: 83418-1-PBS. Standard: Ag34669. Range: 2.5-80 ng/mL.