

For Research Use Only

THOC5 Recombinant antibody

Catalog Number:83435-1-RR



Basic Information

Catalog Number:

83435-1-RR

Size:

800 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6664

GenBank Accession Number:

BC003615

GeneID (NCBI):

8563

UNIPROT ID:

Q13769

Full Name:

THO complex 5

Calculated MW:

79 kDa

Purification Method:

Protein A purification

CloneNo.:

240013F7

Recommended Dilutions:

IF/ICC 1:200-1:800

Applications

Tested Applications:

IF/ICC, FC (Intra), ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : U2OS cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ U2OS cells were intracellularly stained with 0.25 ug Thoc5 Recombinant Antibody (83435-1-RR, Clone:240013F7) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using THOC5 antibody (83435-1-RR, Clone: 240013F7) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).

Immunofluorescent analysis of (4% PFA) fixed U2OS cells using THOC5 antibody (83435-1-RR, Clone: 240013F7) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).