

For Research Use Only

Mouse CD200 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83437-3-PBS

Basic Information

Catalog Number:

83437-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

17470

UNIPROT ID:

O54901

Full Name:

CD200 antigen

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

240391B1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00439-1, MOUSE CD200 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83437-1-PBS. Detection antibody: 83437-3-PBS. Standard: Eg0920. Range: 0.156-20 ng/mL.