

For Research Use Only

TREM2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83438-2-PBS

Basic Information

Catalog Number: 83438-2-PBS	GenBank Accession Number: BC032362	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54209	CloneNo.: 240384D4
Source: Rabbit	UNIPROT ID: Q9NZC2	
Isotype: IgG	Full Name: triggering receptor expressed on myeloid cells 2	
	Calculated MW: 222 aa, 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00440-1, TREM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83438-1-PBS. Detection antibody: 83438-2-PBS. Standard: Eg0747. Range: 0.313-40 ng/mL