

For Research Use Only

TREM2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83438-3-PBS

Basic Information

Catalog Number:

83438-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC032362

GeneID (NCBI):

54209

UNIPROT ID:

Q9NZC2

Full Name:

triggering receptor expressed on myeloid cells 2

Calculated MW:

222 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

240384B12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00440-2, TREM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83438-3-PBS. Detection antibody: 83438-5-PBS. Standard: Eg0747. Range: 0.313-40 ng/mL