

For Research Use Only

# TREM2 Recombinant antibody, PBS Only



Catalog Number: 83438-6-PBS

## Basic Information

Catalog Number:

83438-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC032362

GeneID (NCBI):

54209

UNIPROT ID:

Q9NZC2

Full Name:

triggering receptor expressed on myeloid cells 2

Calculated MW:

222 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

240384B7

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Non-treated and PMA treated THP-1 cells were subjected to SDS PAGE followed by western blot with 83438-6-RR (TREM2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83438-6-PBS in a different storage buffer formulation.