

For Research Use Only

POLR2H Recombinant antibody, PBS Only (Capture)



Catalog Number:83439-2-PBS

Basic Information

Catalog Number: 83439-2-PBS	GenBank Accession Number: BC000739	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5437	CloneNo.: 240410G3
Source: Rabbit	UNIPROT ID: P52434	
Isotype: IgG	Full Name: polymerase (RNA) II (DNA directed) polypeptide H	
Immunogen Catalog Number: AG7166	Calculated MW: 17 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00447-1, POLR2H Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83439-2-PBS. Detection antibody: 83439-1-PBS. Standard: Ag7166. Range: 0.78-100 ng/mL.