

For Research Use Only

SLC2A10 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83441-2-PBS

Basic Information

Catalog Number:

83441-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14210

GenBank Accession Number:

BC101657

GeneID (NCBI):

81031

UNIPROT ID:

O95528

Full Name:

solute carrier family 2 (facilitated glucose transporter), member 10

Calculated MW:

541 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

240335H8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00450-1, SLC2A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83441-2-PBS. Detection antibody: 83441-1-PBS. Standard: Ag14210. Range: 0.78-100 ng/mL.