

For Research Use Only

AGA Recombinant antibody, PBS Only (Capture)



Catalog Number: 83442-5-PBS

Basic Information

Catalog Number: 83442-5-PBS	GenBank Accession Number: BC012392	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 175	CloneNo.: 240396G9
Source: Rabbit	UNIPROT ID: P20933	
Isotype: IgG	Full Name: aspartylglucosaminidase	
Immunogen Catalog Number: AG10681	Calculated MW: 346 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00451-2, AGA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83442-5-PBS. Detection antibody: 83442-3-PBS. Standard: Ag10681. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00451-3, AGA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83442-5-PBS. Detection antibody: 83442-2-PBS. Standard: Ag10681. Range: 0.625-80 ng/mL