

For Research Use Only

NEDD1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83443-1-PBS

Basic Information

Catalog Number:

83443-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33603

GenBank Accession Number:

BC027605

GeneID (NCBI):

121441

UNIPROT ID:

Q8NHV4

Full Name:

neural precursor cell expressed,
developmentally down-regulated 1

Calculated MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

240404B1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00442-2, NEDD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83443-5-PBS. Detection antibody: 83443-1-PBS. Standard: Ag33603. Range: 3.13-100 ng/mL.