

For Research Use Only

NEDD1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83443-5-PBS

Basic Information

Catalog Number: 83443-5-PBS	GenBank Accession Number: BC027605	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 121441	CloneNo.: 240404D3
Source: Rabbit	UNIPROT ID: Q8NHV4	
Isotype: IgG	Full Name: neural precursor cell expressed, developmentally down-regulated 1	
Immunogen Catalog Number: AG33603	Calculated MW: 72 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00442-1, NEDD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83443-5-PBS. Detection antibody: 83443-4-PBS. Standard: Ag33603. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00442-2, NEDD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83443-5-PBS. Detection antibody: 83443-1-PBS. Standard: Ag33603. Range: 3.13-100 ng/mL