

For Research Use Only

SERINC5 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83444-4-PBS

Basic Information

Catalog Number: 83444-4-PBS	GenBank Accession Number: BC101280	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 256987	CloneNo.: 240407H5
Source: Rabbit	UNIPROT ID: Q86VE9	
Isotype: IgG	Full Name: serine incorporator 5	
Immunogen Catalog Number: AG25265	Calculated MW: 423 aa, 47 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00453-2, SERINC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83444-4-PBS. Detection antibody: 83444-3-PBS. Standard: Ag25265. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00453-3, SERINC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83444-1-PBS. Detection antibody: 83444-4-PBS. Standard: Ag25265. Range: 3.13-100 ng/mL