

For Research Use Only

ESR2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83447-1-PBS

Basic Information

Catalog Number: 83447-1-PBS	GenBank Accession Number: BC024181	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2100	CloneNo.: 240412A12
Source: Rabbit	UNIPROT ID: Q92731	
Isotype: IgG	Full Name: estrogen receptor 2 (ER beta)	
Immunogen Catalog Number: AG5103	Calculated MW: 59 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00441-2, ESR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83447-3-PBS. Detection antibody: 83447-1-PBS. Standard: Ag5103. Range: 0.781-50 ng/mL