

For Research Use Only

TMEM119 Recombinant antibody, PBS Only (Capture)



Catalog Number:83448-5-PBS

Basic Information

Catalog Number:

83448-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_181724.2

GeneID (NCBI):

338773

UNIPROT ID:

Q4V9L6

Full Name:

transmembrane protein 119

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

240433F1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00444-1, TMEM119 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83448-5-PBS. Detection antibody: 83448-3-PBS. Standard: Eg0796. Range: 0.313-40 ng/mL.