

For Research Use Only

PLXDC2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83449-2-PBS

Basic Information

Catalog Number: 83449-2-PBS	GenBank Accession Number: BC012885	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84898	CloneNo.: 240174F1
Source: Rabbit	UNIPROT ID: Q6UX71	
Isotype: IgG	Full Name: plexin domain containing 2	
Immunogen Catalog Number: AG2940	Calculated MW: 480 aa, 54 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00445-1, PLXDC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83449-3-PBS. Detection antibody: 83449-2-PBS. Standard: Ag2940. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00445-3, PLXDC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83449-2-PBS. Detection antibody: 83449-1-PBS. Standard: Ag2940. Range: 0.313-40 ng/mL