

For Research Use Only

ABCG8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83450-2-PBS

Basic Information

Catalog Number: 83450-2-PBS	GenBank Accession Number: BC113657	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64241	CloneNo.: 240293E4
Source: Rabbit	UNIPROT ID: Q9H221	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family G (WHITE), member 8	
Immunogen Catalog Number: AG18526	Calculated MW: 673 aa, 76 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00446-2, ABCG8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83450-2-PBS. Detection antibody: 83450-4-PBS. Standard: Ag18526. Range: 1.25-80 ng/mL