

For Research Use Only

# ABCG8 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83450-4-PBS

## Basic Information

Catalog Number:

83450-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18526

GenBank Accession Number:

BC113657

GeneID (NCBI):

64241

UNIPROT ID:

Q9H221

Full Name:

ATP-binding cassette, sub-family G (WHITE), member 8

Calculated MW:

673 aa, 76 kDa

Purification Method:

Protein A purification

CloneNo.:

240293H3

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00446-2, ABCG8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83450-2-PBS. Detection antibody: 83450-4-PBS. Standard: Ag18526. Range: 1.25-80 ng/mL