

For Research Use Only

GNB3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83452-1-PBS

Basic Information

Catalog Number: 83452-1-PBS	GenBank Accession Number: BC000115	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2784	CloneNo.: 240281B5
Source: Rabbit	UNIPROT ID: P16520	
Isotype: IgG	Full Name: guanine nucleotide binding protein (G protein), beta polypeptide 3	
Immunogen Catalog Number: AG0127	Calculated MW: 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00443-2, GNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83452-3-PBS. Detection antibody: 83452-1-PBS. Standard: Ag0127. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00443-3, GNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83452-2-PBS. Detection antibody: 83452-1-PBS. Standard: Ag0127. Range: 0.78-100 ng/mL