

For Research Use Only

HGF Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83453-2-PBS

Basic Information

Catalog Number: 83453-2-PBS	GenBank Accession Number: BC130286	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3082	CloneNo.: 240383E11
Source: Rabbit	UNIPROT ID: P14210	
Isotype: IgG	Full Name: hepatocyte growth factor (hepapoietin A; scatter factor)	
Immunogen Catalog Number: AG25478		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00454-1, HGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83453-2-PBS. Detection antibody: 83453-1-PBS. Standard: Ag25478. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00454-2, HGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83453-3-PBS. Detection antibody: 83453-2-PBS. Standard: Ag25478. Range: 0.313-40 ng/mL