

For Research Use Only

DFNA5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83454-3-PBS

Basic Information

Catalog Number: 83454-3-PBS	GenBank Accession Number: BC019689	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1687	CloneNo.: 240394F1
Source: Rabbit	UNIPROT ID: O60443	
Isotype: IgG	Full Name: deafness, autosomal dominant 5	
Immunogen Catalog Number: AG3746	Calculated MW: 496 aa, 55 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00455-3, DFNA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83454-2-PBS. Detection antibody: 83454-3-PBS. Standard: Ag3746. Range: 0.78-100 ng/mL