

For Research Use Only

IFITM2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83455-3-PBS

Basic Information

Catalog Number: 83455-3-PBS	GenBank Accession Number: BC009696	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10581	CloneNo.: 240318C1
Source: Rabbit	UNIPROT ID: Q01629	
Isotype: IgG	Full Name: interferon induced transmembrane protein 2 (1-8D)	
Immunogen Catalog Number: AG3451	Calculated MW: 132 aa, 15 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00456-2, IFITM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83455-4-PBS. Detection antibody: 83455-3-PBS. Standard: Ag3451. Range: 0.78-100 ng/mL