

For Research Use Only

IFITM2 Recombinant antibody, PBS Only



Catalog Number: 83455-5-PBS

Basic Information

Catalog Number: 83455-5-PBS	GenBank Accession Number: BC009696	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10581	CloneNo.: 240318B6
Source: Rabbit	UNIPROT ID: Q01629	
Isotype: IgG	Full Name: interferon induced transmembrane protein 2 (1-8D)	
Immunogen Catalog Number: AG3451	Calculated MW: 132 aa, 15 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83455-5-RR (IFITM2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83455-5-PBS in a different storage buffer formulation.



1×10^6 MCF-7 cells were intracellularly stained with 0.25 ug IFITM2 Recombinant antibody (83455-5-RR, Clone:240318B6) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 83455-5-PBS in a different storage buffer formulation.