

For Research Use Only

PMS1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83456-4-PBS

Basic Information

Catalog Number:

83456-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28132

GenBank Accession Number:

BC096331

GeneID (NCBI):

5378

UNIPROT ID:

P54277

Full Name:

PMS1 postmeiotic segregation increased 1 (*S. cerevisiae*)

Purification Method:

Protein A purification

CloneNo.:

240408H3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00457-1, PMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83456-3-PBS. Detection antibody: 83456-4-PBS. Standard: Ag28132. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00457-2, PMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83456-2-PBS. Detection antibody: 83456-4-PBS. Standard: Ag28132. Range: 0.156-20 ng/mL