

For Research Use Only

PMS1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83456-7-PBS



Basic Information

Catalog Number: 83456-7-PBS	GenBank Accession Number: BC096331	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5378	CloneNo.: 240408B6
Source: Rabbit	UNIPROT ID: P54277	
Isotype: IgG	Full Name: PMS1 postmeiotic segregation increased 1 (S. cerevisiae)	
Immunogen Catalog Number: AG28132		

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00457-4, Human PMS1 Recombinant Matched Antibody Pair - PBS only. 83456-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag28132. 83456-7-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL

Jurkat cell lysate was measured. The human PMS1 concentration of detected samples was determined to be 6.03 ng/mL (based on a 1.0 mg/mL extract load).