

For Research Use Only

IL1RL1 Recombinant antibody, PBS Only



Catalog Number: 83459-1-PBS

Basic Information

Catalog Number: 83459-1-PBS	GenBank Accession Number: BC030975	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9173	CloneNo.: 240142G3
Source: Rabbit	UNIPROT ID: Q01638-2	
Isotype: IgG	Full Name: interleukin 1 receptor-like 1	
	Calculated MW: 63 kDa, 37 kDa, 30 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data