

For Research Use Only

RNF133 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83462-4-PBS

Basic Information

Catalog Number: 83462-4-PBS	GenBank Accession Number: BC022038	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 168433	CloneNo.: 240397H1
Source: Rabbit	UNIPROT ID: Q8WVZ7	
Isotype: IgG	Full Name: ring finger protein 133	
Immunogen Catalog Number: AG14537	Calculated MW: 376 aa, 42 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00458-1, RNF133 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83462-2-PBS. Detection antibody: 83462-4-PBS. Standard: Ag14537. Range: 0.313-40 ng/mL