

For Research Use Only

TRMT2A Recombinant antibody, PBS Only (Capture)



Catalog Number:83463-1-PBS

Basic Information

Catalog Number:

83463-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9112

GenBank Accession Number:

BC013352

GeneID (NCBI):

27037

UNIPROT ID:

Q8IZ69

Full Name:

TRM2 tRNA methyltransferase 2 homolog A (*S. cerevisiae*)

Calculated MW:

625 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

240444E7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00459-2, TRMT2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83463-1-PBS. Detection antibody: 83463-3-PBS. Standard: Ag9112. Range: 0.625-80 ng/mL.