

For Research Use Only

TRMT2A Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83463-5-PBS

Basic Information

Catalog Number: 83463-5-PBS	GenBank Accession Number: BC013352	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 27037	CloneNo.: 240444F6
Source: Rabbit	UNIPROT ID: Q8IZ69	
Isotype: IgG	Full Name: TRM2 tRNA methyltransferase 2 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG9112	Calculated MW: 625 aa, 69 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00459-1, TRMT2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83463-5-PBS. Detection antibody: 83463-4-PBS. Standard: Ag9112. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00459-3, TRMT2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83463-2-PBS. Detection antibody: 83463-5-PBS. Standard: Ag9112. Range: 0.625-80 ng/mL.