

For Research Use Only

LSM14A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83464-5-PBS

Basic Information

Catalog Number: 83464-5-PBS	GenBank Accession Number: BC016842	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26065	CloneNo.: 240380D2
Source: Rabbit	UNIPROT ID: Q8ND56	
Isotype: IgG	Full Name: LSM14A, SCD6 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG12997	Calculated MW: 463 aa, 51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00463-1, LSM14A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83464-4-PBS. Detection antibody: 83464-5-PBS. Standard: Ag12997. Range: 0.625-80 ng/mL.