

For Research Use Only

GRIN2A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83465-1-PBS

Basic Information

Catalog Number: 83465-1-PBS	GenBank Accession Number: NM_000833	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2903	CloneNo.: 240445B5
Source: Rabbit	UNIPROT ID: Q12879	
Isotype: IgG	Full Name: glutamate receptor, ionotropic, N-methyl D-aspartate 2A	
Immunogen Catalog Number: AG29101	Calculated MW: 165 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00461-3, GRIN2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83465-4-PBS. Detection antibody: 83465-1-PBS. Standard: Ag29101. Range: 0.313-40 ng/mL