

For Research Use Only

HSPA2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83466-3-PBS

Basic Information

Catalog Number: 83466-3-PBS	GenBank Accession Number: BC036107	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3306	CloneNo.: 240421A4
Source: Rabbit	UNIPROT ID: P54652	
Isotype: IgG	Full Name: heat shock 70kDa protein 2	
Immunogen Catalog Number: AG3539	Calculated MW: 639 aa, 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00465-1, HSPA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83466-1-PBS. Detection antibody: 83466-3-PBS. Standard: Ag3539. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00465-3, HSPA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83466-3-PBS. Detection antibody: 83466-2-PBS. Standard: Ag3539. Range: 1.25-80 ng/mL