

For Research Use Only

ELAC2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83468-1-PBS

Basic Information

Catalog Number: 83468-1-PBS	GenBank Accession Number: BC004158	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 60528	CloneNo.: 240487F4
Source: Rabbit	UNIPROT ID: Q9BQ52	
Isotype: IgG	Full Name: elaC homolog 2 (E. coli)	
Immunogen Catalog Number: AG0105	Calculated MW: 92 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00471-1, ELAC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83468-4-PBS. Detection antibody: 83468-1-PBS. Standard: Ag0105. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00471-2, ELAC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83468-2-PBS. Detection antibody: 83468-1-PBS. Standard: Ag0105. Range: 0.313-40 ng/mL